nuvitamed

- IT** Termometro digitale frontale / auricolare
- ES** Termómetro digital para frente/oído
- EN** Digital forehead / ear thermometer
- FR** Thermomètre frontal/auriculaire numérique
- DE** Digitales Stirn-/Ohrthermometer

2075

- IT** Manuale di istruzioni
- ES** Manual de instrucciones de uso
- EN** Instruction manual
- FR** Manuel d'instructions d'utilisation
- DE** Bedienungsanleitung

AOJ-20A

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Thank you for purchasing this Forehead & Ear Infrared Digital Thermometer. Please read the Operating Manual carefully before using the thermometer, and put it in a safe and secure place for reference.

The infrared thermometer can be used to read the body temperature by measuring the ear and forehead temperature, which is suitable for adult and children (ear test mode only for the child above 3 months).

WARNINGS AND PRECAUTIONS

1. Keep out of reach of children under 12 years.
2. Never immerse the thermometer into water or other liquids (not waterproof). For cleaning and disinfecting, please follow the instructions in the "Care and cleaning" section.
3. Never use the thermometer for purposes other than those it has been intended for. Please follow the general safety precautions when using on children.
4. Keep the thermometer away from direct exposure to the sun and keep it in a dust-free, dry area, well-ventilated place at a temperature between 5°C (41 °F) - 40°C (104°F) . Do not use the thermometer in high humidity environments. (>95% RH)
5. Do not use the thermometer if there are signs of damage on the measuring sensor or on the instrument itself. If damaged, do not attempt to repair the instrument! Please contact the dealer.
6. This thermometer consists of high-quality precision parts. Do not drop the instrument. Protect it from severe impact and shock. Do not twist the instrument or the measuring sensor. Please consult your doctor if you see symptoms such as unexplained irritability, vomiting, diarrhea, dehydration, changes in appetite or activity, seizure, muscle pain, shivering, stiff neck, pain when urinating, etc., even in the absence of fever.
7. Even in the absence of fever, those who exhibit a normal temperature may still need to receive medical attention. People who are on antibiotics, analgesics, or antipyretics should not be assessed solely on temperature readings to determine the severity of their illness.
8. Temperature elevation may signal a serious illness, especially in adults who are old, frail, have a weakened immune system, or neonates and infants. Please seek professional advice immediately when there is a temperature elevation and if you are taking temperature for whom are:
 - Over 60 years of age (Fever may be blunted or even absent in elderly patients)
 - Having diabetes mellitus or a weakened immune system (e.g., HIV positive, cancer, chemotherapy, chronic steroid treatment, splenectomy)
 - Bedridden (e.g., nursing home patient, stroke, chronic illness)
 - A transplant patient (e.g., liver, heart, lung, kidney)
9. This thermometer is not intended for preterm babies or small-for-gestational age

- babies. This thermometer is not intended to interpret hypothermic temperatures.
10. Do not allow children to take their temperatures unattended.
 11. Use of this thermometer is not intended as a substitute for consultation with your physician or pediatrician.
 12. Clean the thermometer probe after each use.
 13. Do not use the thermometer on newborns or for continuous temperature monitoring purposes.
 14. Do not take a measurement while or immediately after nursing a baby.
 15. Patients should not drink, eat, or be physically active before/while taking the measurement.

PRODUCT DESCRIPTION

OVERVIEW

Infrared Thermometer measures the body temperature based on the infrared energy emitted from the eardrum or the forehead. Users can quickly get measurement results after positioning properly the temperature probe in the ear canal or forehead. Normal body temperature is a range. The following tables shows that this normal range also varies by site. Therefore, readings from different site should not be directly compared. Tell your doctor what type of thermometer you used to take your temperature and on what part of the body. Also bear this in mind if you are diagnosing yourself.

Measuring part	Normal temperature range
Forehead temperature	36.1°C - 37.5°C (97°F - 99.5°F)
Ear temperature	35.8°C - 38°C (96.4°F - 100.4°F)
Oral temperature	35.5°C - 37.5°C (95.9°F - 99.5°F)
Rectal temperature	36.6°C - 38°C (97.9°F - 100.4°F)
Axillary temperature	34.7°C - 37.3°C (94.5°F - 99.1°F)

STRUCTURE

The device consists of enclosure, button, a temperature sensor, Infrared temperature-measuring element, microcomputer-controlled circuit, LCD, backlight and a buzzer.

OPERATING PRINCIPLE

In the nature, for all objects whose temperature is higher than the absolute zero (-273.15°C), for the sake of molecule's thermal motion, radiate electromagnetic wave to the surrounding ambient including the infrared wave without a break, the relationship between the density of radiation energy with the object's temperature conform to the radiation law. The infrared thermometer's working principle is based on the fourth-power law: infer objects' radiation temperature by measuring objects' radiation infrared energy. Thermopile sensor can convert the infrared energy into thermoelectricity, and output as a detected signal after signal processing. The infrared temperature sensor



Figure 1

1. Subject
2. Object lens
3. Heat plate
4. Thermocouple
5. Eye Lens

detects infrared energy emitted by the eardrum. A built-in lens focuses the collected energy, which is then converted into a temperature reading by the thermopiles and measurement circuits.

A. Infrared thermometer probe; **B.** Infrared ray; **C.** Ear canal; **D.** Eardrum; **E.** Computer; **F.** Display; **G.** IR Sensor; **H.** Optical component; **I.** Ear canal; **J.** Tympanic membrane; **K.** External ear; **L.** Probe cover; **M.** Reference.

INDICATIONS FOR USE

The infrared thermometer takes human body temperature via the eardrum or forehead. It applies to all age groups except for babies under three months. Both devices apply to both professional use and home use.

CONTRAINDICATIONS


Do not use the thermometer if the ear is infected with otitis or suppuration.

FEATURES

- Magnetic probe cover and automatic mode switch between forehead & ear mode (OPTIONAL)
- Quick measurement, less than 1 second
- Small Body with super large font
- Multi-functional, can measure ear, forehead
- Easy operation, one button design, to measure ear, forehead
- 40 sets of memories, easy to recall. Switching between mute and un-mute mode
- Audio and visual fever warning. Switching between °C and °F
- Auto shut-down and power-saving
- Backlight color temperature indication
- Patent technology to realize separate measurement of adult and child


PRODUCT STRUCTURE

1. LCD display screen
2. Mode/Memory button
3. Measure/Mute-unmute button
4. Probe
5. Battery cover
6. Probe cover (take it off when measuring ear temperature)

DISPLAY DESCRIPTION

1. Adult forehead mode
2. Child forehead mode (under 12 years old)
3. Ear mode
4. Adjust mode
5. Fahrenheit / Celsius degrees
6. Previous temperature value
7. Low battery indication
8. Mute /Un-mute icon
9. Current temperature value

HOW TO USE YOUR THERMOMETER

IN ORDER TO AVOID THE INACCURACY:

1. Please ensure that the device is used only in the room, and there is no strong wind presence.
2. Please make sure that there is no intense emotion and movement before measuring.
3. If the device is transferred from one condition to another, which has different ambient temperature, it is suggested to deposit for more than 30 minutes.
4. If the tester is transferred from one condition to another, which has different ambient temperature, it is suggested to have a rest for more than 10 minutes.
5. Do not hold the device for a long time, as it is highly sensitive to heat. The device has undergone clinical test, it is safe and accurate when using in accordance with operation manual.

TAKE YOUR FOREHEAD TEMPERATURE

When using the thermometer for the first time, please load the batteries.

Press and release the Measure/Mute-unmute button in one second to power on. Press the Mode/ Memory button to choose "Adult forehead mode" or "Child forehead mode". Then point the temperature probe at the center of the forehead, with a distance of 0-3 cm (0-1.18 inches).

Press and release the Measure/Mute-unmute button in 1 second. The beep is heard, you can now read the value.

NOTE: The forehead measurement is an indicative reading. The measured forehead temperature can fluctuate up to 1 °F/0.5 °C from your actual body temperature. Please be aware of the factors that influence the accuracy as described in the section "Temperature taking tips" and "WARNINGS AND PRECAUTIONS".

- ⚠ If the eyebrow area is covered with hair, sweat or dirt, please clean the area beforehand to improve the reading accuracy.
- ⚠ Always check if the lens is clean.
- ⚠ Always make sure the user and the thermometer will have been in the same room for at least 30 minutes prior to the measurement.

TAKE YOUR EAR TEMPERATURE

1. Press the Measure/Mute-unmute button to power on.
2. Press the Mode/Memory button to choose "Ear mode" and take off the probe cover (if with magnetic probe cover, it is automatically switched between forehead & ear mode), fit the probe snugly into the ear canal.
3. Press and release the Measure/Mute-unmute button in 1 second, the beep is heard, you can now read the value.

- ⚠ Please make sure that the ear is clean, with no earwax or obstructions.
- ⚠ The right ear reading may differ from the reading taken at the left ear. Therefore, always take the temperature in the same ear.


NOTE:

Children under 1 year: Pull the ear straight back.

Children aged 1 year to adult: Pull the ear up and back.

- ⚠ Do not force the thermometer into the ear canal. Otherwise, the ear canal may get injured.
- ⚠ When taking the temperature on an adult, gently pull the ear up and back to make sure the ear canal is straight, so that the temperature probe can receive an infrared ray from the eardrum.
- ⚠ Be careful when taking temperature on a child, whose ear canal is small.

AFTER A MEASUREMENT

-  Once the reading has been completed, remove the thermometer away from the forehead/ear and observe temperature. After each measurement, you can enter the recall mode and query earlier temperature readings. Do not hold the thermometer for a long time, because it is sensitive to the ambient temperature.
-  After each measurement, clean the temperature probe with a soft cloth, and put the thermometer in a dry and well-ventilated place. You should wait at least 10 seconds between each measurement.
-  It is dangerous to make a self-diagnosis or self-treatment based on the obtained measurement results. For such purposes, please consult a doctor.

READ YOUR TEMPERATURE

"T" indicates a temperature reading.

In forehead or ear mode. If $32^{\circ}\text{C} \leq T \leq 37.5^{\circ}\text{C}$ ($89.6^{\circ}\text{F} \leq T \leq 99.5^{\circ}\text{F}$), the green light will be displayed, with one long beep. If $37.6^{\circ}\text{C} \leq T \leq 38.0^{\circ}\text{C}$ ($99.7^{\circ}\text{F} \leq T \leq 100.4^{\circ}\text{F}$), the orange light will be displayed with 6 short beeps, which is a warning that you may have a low fever.

If $38.1^{\circ}\text{C} \leq T \leq 42.9^{\circ}\text{C}$ ($100.6^{\circ}\text{F} \leq T \leq 109.2^{\circ}\text{F}$), the red light will be displayed with 6 short beeps, which is a warning that you may have a high fever.

In Adjust mode. If $0^{\circ}\text{C} \leq T \leq 100^{\circ}\text{C}$ ($32.0^{\circ}\text{F} \leq T \leq 199.0^{\circ}\text{F}$), the white light will be displayed with one long beep. This mode is not for human temperature measurement, for maintenance use only.

SWITCHING BETWEEN MUTE AND UN-MUTE

When the thermometer is turned on, keep pressing the Measure/Mute-unmute for 5 seconds, to switch from un-mute to mute.

CHECKING 40 SETS OF MEMORY DATA

When the thermometer is turned on or off, keep pressing the Mode/ Memory button for 5 seconds to go to the memory mode, press this button again to check the 40 sets of memories one by one. If no value, it will display "---".

°C/°F CONVERSION

Open the battery cover, use the toggle switch to change the °C/°F.

AUTOMATICALLY TURN OFF

The unit will shut down automatically after 13 seconds of no use.

CAUTION

All memory records will be lost when uninstall or reinstall the battery.

REPLACE THE BATTERY

-  Put two AAA batteries correctly into the compartment.
- Remove the batteries if the thermometer will not be used for more than two months.

TEMPERATURE TAKING TIPS

1. It is important to know each individual's normal temperature when they are well. This is the only way to accurately diagnose a fever. Record readings twice a day (early morning and late afternoon). Take the average of the two temperatures to calculate normal oral equivalent temperature. Always take the temperature in the same location, since the temperature readings may vary from different locations on the forehead.
2. A child's normal temperature can be as high as 99.9°F (37.7°C) or as low as 97.0°F (36.1°C). Please note that this unit reads 0.5°C (0.9°F) lower than a rectal digital thermometer.
3. External factors may influence ear temperatures, including when an individual has:
 - been lying on one ear or the other
 - had their ears covered
 - been exposed to very hot or very cold temperatures
 - been recently swimming or bathing
4. In these cases, remove the individual from the situation and wait 20 minutes prior to taking a temperature. Use the untreated ear if prescription ear drops or other ear medications have been placed in the ear canal.
5. Holding the thermometer for too long in the hand before taking a measurement can cause the device to warm up. This means the measurement could be incorrect.
6. Patients and the thermometer should stay in steady-state room condition for at least 30 minutes.
7. Before placing the thermometer sensor onto the forehead, remove dirt, hair, or sweat from the forehead area. Wait 10 minutes after cleaning before taking measurement.
8. Use an alcohol swab to carefully clean the sensor and wait for 5 minutes before taking a measurement on another patient. Wiping the forehead with a warm or cool cloth may impact your reading. It is advised to wait 10 minutes before taking a reading.
9. In the following situations it is recommended that 3-5 measurements in the same location be taken and the highest one taken as the reading:
 - Newborn infants in the first 100 days.
 - Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 - When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the instrument and obtains consistent readings.

CARE AND CLEANING

Use an alcohol swab or cotton swab moistened with 70% alcohol to clean the thermometer casing and the measuring probe. After the alcohol has completely dried out, you can take a new measurement. Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the instrument in water or other cleaning liquids. Take care not to scratch the surface of the LCD screen.

ERROR AND TROUBLESHOOTING

Symptom	Possible cause	Description & solution
Failed to power on	The battery level is too low.	Replace with a new battery
	Polarities of the batteries are reversed.	Ensure the batteries are in the right position
	The thermometer is damaged	Contact dealer
The reading is too low	The lens of the probe is dirty.	Clean the lens with a cotton swab.
	The distance of the item and target is too far	Keep in contact with forehead, or put the probe into the Ear Canal.
	You have just come from a cold environment	Stay in a warmer room for at least 30 minutes before taking a reading
The reading is too high	You have just come from a hot environment.	Stay in an adequately cool room for at least 30 minutes before taking a reading
	The ambient temperature is not in range.	One long beep, and red backlit for 3 seconds. Take a measurement under an ambient temperature between 5°C (41.0°F) and 40 °C (104°F).
	Memory error or calibration is not finished. Item will be turned off automatically after 3 seconds.	One long beep, and red backlit for 3 seconds. Contact dealer.
	In ear or forehead mode, $T > 42.9^{\circ}\text{C}$ (109.2°F). In adjust mode, $T > 100^{\circ}\text{C}$ (199°F)	One long beep, the screen will be displayed in red.

	<p>In ear or forehead mode, $T < 32^{\circ}\text{C}$ (89.6°F). In adjust mode, $T < 0^{\circ}\text{C}$ (32°F)</p>	<p>One long beep, the screen will be displayed in red.</p>
	<p>The power voltage is between 2.4V-2.7V</p>	<p>The battery icon will flicker, it means the battery is low, but you still can use the thermometer.</p>
	<p>The power voltage is lower than 2.4V</p>	<p>The battery icon will flicker, the item will be turned off automatically after 3seconds. Please replace with a new battery</p>

TECHNICAL SPECIFICATIONS

Product name	Forehead & Ear in fared thermometer
Applicable regulations and laws	ASTM E 1965-98
	ISO80601-2-56
Display	Segment LCD, Four colour LED backlight (white, green, orange, red)
Temperature units	$^{\circ}\text{C}$ / $^{\circ}\text{F}$, switchable
Power supply	2x1.5V AAA Batteries
Measurement range	Measurement: 32.0°C – 42.9°C (89.6°F – 109.2°F)
	Direct mode: 0°C – 100°C (32°F – 199°F)
Accuracy (Laboratory)	Measurement mode $\pm 0.2^{\circ}\text{C}$ / $\pm 0.4^{\circ}\text{F}$
	Direct mode $\pm 1.0^{\circ}\text{C}$ / 2.0°F
Display resolution	0.1°C / $^{\circ}\text{F}$
Shutdown Memory	40 groups of measured temperature.
Operational conditions	Temperature: 5°C – 40°C (41°F – 104°F). Humidity: 15-95%RH, non-condensing. Atmospheric pressure: 70-106kPa.

Storage condition	Ambient Temperature: -20°C-55°C(-4°F-131°F). Relative Humidity: 0-95% RH, non-condensing. Atmospheric pressure: 50kPa to 106kPa.
Battery	2*AAA, can be used for more than 3000 times
Weight & Dimension	60g (without battery), 143x35x41mm

SYMBOLS

Symbol	Description
	Type BF applied part.
	Information about a manufacturer.
	Please read the instructions carefully.
	Waste electrical materials should be sent to a dedicated collection point for recycling.
	Production Date
LOT	Batch number
	IMPORTANT Inaccurate reading or thermometer damage may occur if the thermometer is not correctly used.
IP22	2 Protected against solid foreign objects of 12,5 mm Ø and greater; 2 If keep the thermometer in 15 degree angle, it still can prevent the water drop.
	CE marking
	Authorized Representative In The European Community
	Keep away from sunlight
	Fragile, handle with care
	Keep dry
	This way up

APPENDIX: EMC INFORMATION-GUIDANCE AND MANUFACTURER'S

DECLARATIONS CAUTION: • The Infrared Thermometer needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided for in the ACCOMPANYING DOCUMENTS. • Portable and mobile RF communications equipment can affect Infrared Thermometer • The Infrared Thermometer should not be used adjacent to or stacked with other equipment. • Guidance and manufacturer's declaration – electromagnetic emission –for all EQUIPMENT AND SYSTEMS.

Guidance and manufacturer's statement- Electromagnetic emission	
The Infrared Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared Thermometer should assure that it is used in such an environment.	
Emission test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions IEC61000-3-3	Not applicable

Guidance and manufacturer's declaration -Electromagnetic Immunity		
The Infrared Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared Thermometer should assure that it is used in such an environment.		
Immunity Test	IEC 60601-1-2 Test Level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 KV contact ±2 Kv, ±4 kV, ±8 kV, ±15 kV air	±8 KV contact ±2 Kv, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable
Power frequency magnetic field IEC 61000-4-8	30 A/m 50HZ / 60Hz	30 A/m 50HZ / 60Hz
Conducted RF IEC61000-4-6	Not applicable	Not applicable
Radiated RF IEC61000-4-3	10 V/m; 80MHz - 2,7 GHz 80 % AM at 1 kHz	10 V/m; 80MHz - 2,7 GHz 80 % AM at 1 kHz
NOTE UT is the a.c. mains voltage prior to application of the test level		

ELECTROMAGNETIC IMMUNITY

Guidance and manufacturer's declaration - Electromagnetic Immunity

The Infrared Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared Thermometer should assure that it is used in such an environment.

Test Frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
385	380-390	TETRA	Pulse modulation 18 Hz	1.8	0.3	27	27
400	430-470	GMR5-460 FRS 460	FM±5 kHz deviation 1 kHz sine	2	0.3	28	28
710							
745	704-787	LTE Band 13,17	Pulse modulation 217 Hz	0.2	0.3	9	9
780							
810							
870	800-960	GSM 800/900, TETRA 400, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28	28
930							
810							
870	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28	28
930							
2450	2400-2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450,	Pulse modulation 217 Hz	2	0.3	28	28
2450	2400-2570	Bluetooth WLAN 802.11 b/g/n, RFID 2450,	Pulse modulation 217 Hz	2	0.3	28	28
810							
870	5100-5800	LTE Band 7 WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9	9
930							

Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)

WARRANTY - TERMS AND CONDITIONS

This product benefits of a warranty of 24 months on material and manufacturing defects, starting from the date of purchase (see sales slip). The warranty of 24 months does not include damages caused by the usual use of parties identified as "consumable" (e.g., batteries, brush heads or parts subject to usury). The legal guarantee of 24 months is void if: **1.** The product has undergone aesthetic damage due to improper use not in accordance with instructions in the manual. **2.** This product has been modified and/or tampered with. **3.** The cause of the failure was due to poor maintenance of the individual components and/or accessories and/or supplies (e.g. oxidation and/or scaling due to the retention of water or other liquids, sediment blocking the sensor, leak of corrosive liquid from batteries). The following is excluded from the legal guarantee of 24 months: **1.** Costs related to replacements and/or repair of parts subject to wear or costs for ordinary maintenance of the product. **2.** The costs and risks involved in transporting the product to and from the store where you purchased or otherwise authorized collection centre to receive the products under warranty. **3.** Damage caused by or resulting from improper installation or improper use not in accordance with the directions in the instruction manual. **4.** Damage due to natural disasters, accidental events or adverse conditions not compatible with the product. **5.** Defects that have a negligible effect on product performance. The manufacturer, distributor and all the parties involved in the sale do not assume any liability for losses and economic damage from any malfunction of the product. In accordance with current regulations the manufacturer, distributor and all the parties involved in the sale are not responding in any case for damages, including direct, indirect ones, loss of net income, loss of savings and additional damage and other details consequences going beyond the damage caused by the breach of warranty, contract, strict liability, wrongdoing or due to other causes, resulting from the use or inability to use the product and/or paper and electronic documents, including the lack of service. For further information on the help service visit the website www.nuvitababy.com

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All rights reserved
Share Info GmbH
Shenzhen AOJ Medical Technology Co., Ltd.
Anteprima Brands International Ltd