# nuvitamed

- IT Termometro Digitale Auricolare
- ES Termómetro Digital de Oído
- EN Digital Ear Thermometer
- FR Thermomètre Auriculaire Digital
- DE Digitales Ohrthermometer

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
Congratulations on your purchase of this digital thermometer that allows. For a proper and in compliance with the safety of the user to apply strictly the precautions and warnings listed in this manual. Before using the thermometer read carefully and thoroughly this manual and retain it for future reference.

## PRODUCT DESCRIPTION

1. Ear Probe
2. LED Display
3. ON/Memory button
4. Battery Cover
5. Thermometer
6. SCAN button

1. Digital Readout
2. Cry Indicator-Fever
3. Smile Indicator- Normal
4. Measure indicator
5. Ear Temperature Indicator
6. Unit (°C/°F)
7. Memories Indicator
8. Memories Number
9. Low Battery Indicator

### ATTENTION

Clean the probe with alcohol after every use to maintain hygiene and accuracy.

## PREPARATION

Make sure the batteries are charged. If necessary replace them with two new AAA batteries. To obtain reliable measurements wait 30 minutes to adjust the thermometer at room temperature. The temperature changes can cause unreliable measurements, so do not place the thermometer in front of air conditioners or heat radiators before the measurement because you can get inconsistent results. Make sure that the ear is clean and dry since the presence of ear wax or residues of epithelium can lead to inconsistent results. To ensure accurate results the measurements are made at rest because the body temperature may rise after exercise or after a bath. During measurement, do not talk or move your body. If the measurement is incorrect, you can measure it again. However, take a 5-15 minutes rest before the next measurement.

## HOW TO MEASURE THE TEMPERATURE

To measure the temperature from the ear:

- Press [ON] button to switch on the thermometer.
- Gently insert the Probe into the ear
- Short Press the [SCAN] button for measurement.
- Keep the probe in ear until you hear a beep sound
- LCD screen will show the temperature
- When measuring, hold the thermometer firmly and not touch any button.


The result will be automatically recorded in the memory of the thermometer and assigned the highest memory number. If in the memory of the thermometer there are already 10 records, the oldest is deleted. You should not run more than 3 consecutive measurements. Before performing any 4th measurement please wait 10 minutes. Upon completion and after about 30 seconds of inactivity, the thermometer automatically switches to clock mode. If the sleep function is active, after about 1 minute of inactivity the screen goes off automatically

### MEASUREMENT UNITS

Press and hold the [SCAN] button for 10 seconds, keep holding [SCAN] button and press the [ON] button for 1 second to change the units from Celsius (°C) to Fahrenheit (°F) and vice versa.

### MEMORY

When the thermometer is switched on, press [ON] to see the measurements stored in memory. If you receive a "-" means that there is no data in the memory. The max capacity of memory is 10.

### AUTO-OFF

The thermometer will turn off automatically after no activity for 20 seconds.

### SUGGESTIONS

For non-technical users should be performed 3 measurements test to practice with the device. Consider the maximum value as a continuous use can lead to an error margin of  $\pm 0,3^{\circ}\text{C}$  -  $\pm 0,5^{\circ}\text{F}$ . In healthy individuals there is no difference between left or right ear, so do the tests measure both ears and always measure the temperature from the ear that returns the highest value.

## **BATTERIES REPLACEMENT**

Remove the battery cover. Always use new batteries of the same type (2 x AAA) and check the polarity in the battery compartment (symbols [+] and [-]). Replace the battery cover.

## **MAINTENANCE**

Clean the device with a soft, dry cloth. Do not immerse the device in liquid, solvents or detergents.

## **STORAGE**

Clean the thermometer and cover the sensor with the lid after each use. Do not expose the thermometer to direct sunlight, high temperatures, humidity, fire, flames and prevent vibration or shock. If you do not use the thermometer for a very long period, remove the batteries.

## **CALIBRATION**

This thermometer was calibrated at the factory at the time of construction. If the thermometer is used in accordance with the instructions in this manual is not necessary to recalibrate periodically. In the case of incorrect measurements strongly contact your local dealer for assistance procedure. Do not grope to disassemble, repair or modify the thermometer.

## **CLEANING**



Keep the sensor surface clean otherwise the measurement may be unreliable. Clean the surface of the sensor with a cotton swab moistened with water or alcohol before again using the thermometer to make sure that the sensor surface is dry. If one or more items are missing or damaged, please refer to your local dealer and request the replacement.

## **INTENDED USE**

This digital thermometer is a device for measuring body temperature in domestic use, where continuous monitoring is not required by medical and paramedical staff. The device can be used to measure body temperature from ear and must be used carefully following the instructions in this manual. Any use different from that stated in this manual is considered improper.

## **SAFETY PRECAUTIONS**


The warning symbols and logos appearing in this manual are intended to alert the user to prevent dangerous situations and avoid damage to the user and to others. The following table provides a description:

	Serial number
	Operating temperature
	Operating and storage temperature
	Operating and atmospheric pressure
	This symbol indicates instructions or directions that must be followed in order to avoid damages to user or other people.
	This symbol warns you if you do not comply with these instructions is likely to cause harm to people or objects such as housing, pets or other family members.
	This symbol indicates prohibition of the use of certain objects or accessories specifically mentioned in the instructions.
	This symbol indicates important obligations to be followed explicitly in the instructions.
	This symbol indicates that the product complies with the requirements of protection against electric shock device of type BF.
	Symbol for "MANUFACTURER"
	Symbol for "EUROPEAN REPRESENTATION"
	Consult instructions for use
<b>IP22</b>	IP22: The first number 2: Protected against solid foreign objects of 12,5 mm Ø and greater. The second number: Protected against vertically falling water drops when enclosure tilted up to 15°.
	WARNING: The symbol on this product means that it's an electronic product and following the European directive 2000/96/CE the electronic products have to be dispose on your local recycling centre for safe treatment.
 0598	This device meets the essential requirements and other relevant conditions of the applicable European directives. Complies with MDD 93/42/EEC. 0598 is the Notified Body Number. Model: HTD8208C.

## DISPOSAL

- Used batteries should not be disposed of in the household rubbish. Used Batteries should be deposited at a collection point.
- At the end of its life, the appliance should not be disposed of in household rubbish. Enquire about the options for environment-friendly and appropriate disposal. Take local regulations into account.

## TROUBLESHOOTING

Display	Situation	Solution
	The temperature is higher than 42.9°C (109.2°F)	Make sure the right measuring object
	The temperature is lower than 34°C (93.2°F)	Make sure the right measuring object
	Wrong operating ambient temperature	Make sure the device is in 15-35°C (59~95°F) operating conditions
	System mistake	System mistake happen, the device will turn off after show all the digits and signals.
	Battery capacity is too low. Taking temperature is not allowed.	Please change to new battery
	The operating ambient temperature is not stable	Keep the device ion stable operating ambient temperature more than 30 minutes before measuring.
	Blank display	Make sure the correct batteries are installed and its positive and negative charges.

### WARNINGS

Based on the results obtained by measuring body temperature with this device, it is dangerous for people to make diagnoses and take medicines without consulting your doctor before. A wrong diagnosis and treatment can lead to severe consequences and worsening of the disease. Always follow the advice of your doctor.

Do not touch or swallow the infrared sensor. A dirt sensor can return values of body temperature unreliable and untrustworthy.

In case of body temperatures above 38,0°C and in presence of symptoms of pain refer immediately to your doctor.

Keep away from children because there are small parts that can be swallowed and cause

dangerous situations. The ingestion of any of these parties, such as batteries or the cover sensor, requires immediate intervention by emergency personnel or your doctor. The lack of caution in measuring the temperature through the ear can cause damage and injuries.

This device is not intended to be used by people, including children with physical disabilities, sensory and mental disabilities, or lack of experience or knowledge if not properly taught and trained to use the device by a competent person responsible for their safety.

In the case of illnesses such as ear infections or diseases of the ear drum, use only the thermometer to measure temperature by the forehead. Use the thermometer to measure temperature by the ear could worsen the disease.

If a bath or shower, do not immediately measuring the temperature, since the presence of some water may cause damage to the ear.

If necessary, very gently clean the infrared sensor with a cotton swab moistened with water or alcohol and wait until the surface is completely dry before proceeding with a new use. The use of paper towels can cause scratches on the sensor thereby compromising the efficiency and reliability.

For a correct measurement of body temperature, do not perform the measurement immediately after performing physical activities or after a meal. Wait at least 30 minutes to take the measurement at rest.

The use of methods such as sponging with cold water or ice pack lowers the temperature, so you need not perform the measurement immediately otherwise you get lower values.

Before measuring the body temperature to wait about 30 minutes as the different conditions that we may be one of the rooms, which may affect the results, making the measurement unreliable.

Do not dispose of batteries in a fire as they may explode. Once they no longer charge the batteries disposed of in special containers in accordance with the regulations

The thermometer is not waterproof, so do not immerse in liquid or abrasive cleaners. Do not expose the thermometer to direct sunlight, strong heat or flame.

Once they no longer charge the batteries disposed of in special containers in accordance with the regulations. This digital thermometer is a device for measuring body temperature in domestic use, where continuous monitoring is not required by medical and paramedical staff.

## RECOMMENDATIONS

Do not use the device for purposes other than those indicated in the intended use. Keep the device is affected by falls or impacts. When measuring, do not use a mobile or cordless phone near the device. If unused for a long period, remove the batteries to avoid any leakage of corrosive substances.

## GENERAL INFORMATION

Body temperature does not have a single and absolute value, it can vary for several reasons. The age and sex of the person, for example, can affect health and optimal body temperature can vary during the day going from a minimum of 35,5°C (95,9°F) up to a maximum of 37,8°C (100,0°F). It is recommended to perform measurements under conditions of optimal health at different times of day to get reference values. The digital thermometer detects the infrared heat of the blood flow through the ear canal or through the brow area which, being sprayed can detect temporal artery in a fast, accurate and reliable body temperature. In a healthy person who enjoys good health, the temperature should not vary much between the right and left ear. To find the right one, it is recommended to perform measurements in both ears and always use the one that returns a higher temperature.

## SPECIFICATIONS

---

<b>Operating mode:</b>	Adjusted mode (Body mode)
<b>Measuring mode:</b>	Ear
<b>Reference body site:</b>	Axillary
<b>Product Name:</b>	Nuvita 2071 HTD8208C
<b>Power Supply:</b>	3 VDC (2 x AAA batteries)
<b>Temp. sensing part:</b>	IR sensor
<b>Measuring range:</b>	34,0°C – 42,9°C
<b>Alarm Range:</b>	34.0°C ~ 37.3°C:Green (normal); 37.4°C ~ 38.0°C:Yellow (slight fever); 38.1°C ~ 42.9°C:Red (high fever).
<b>Measuring accuracy:</b>	34.0°C ~ 34.9°C:±0.3°C; 35.0°C ~ 42.0°C:±0.2°C; 42.1°C ~ 42.9°C:±0.3°C;
<b>Dimensions:</b>	160 x 36 x 45 mm
<b>Weight:</b>	60 gr. without batteries
<b>Working environment:</b>	15°C – 35°C
<b>Relative humidity:</b>	20%RH – 80%RH
<b>Storage environment:</b>	-20°C – 55°C
<b>Relative humidity:</b>	20%RH – 90%RH
<b>Atmospheric pressure:</b>	70-106Kpa
<b>Measure time:</b>	<2 seconds
<b>Auto shut off:</b>	≤35 seconds

## COMPLIANCE

Item	Complaint with
Equipment classification	Safety Standards: EN 60601-1: 2006+A1:2013, EN60601-1-2:2015
Type of protection	Internally powered equipment (on battery power)
Degree of protection	Type BF -Applied part(the probe)
Front panel and case labeling	EN ISO 15223-1: 2016
Temperature	EN ISO 80601-2-56: 2012
Home healthcare environment	EN 60601-1-11:2015
Note: Not intended to be sterilized. Not for use in an OXYGEN RICH ENVIRONMENT.	

Accessory name	Number	Item number
AAA battery	2	4200 1511 0600
User manual	1	ZHTF-CE-02-015

Calculated values of the indicators according to EN ISO80601-2-56:2012				
	Group A1	Group A2	Group B	Group C
Bias	0.125	0.180	0.194	0.160
Standard Deviation	0.121	0.130	0.125	0.132
Limits of agreement	0.242	0.259	0.250	0.265
Clinical repeatability (for All Group)	0.097			

Safety classification of ME EQUIPMENT	
Protection against electric shock	Internally powered ME equipment
Degree of protection	Type BF -Applied part (the probe)
Protection against harmful ingress of water or particulate matter	IP22
Mode of operation	Continuous operation
Note: Not intended to be sterilized. Not for use in an OXYGEN RICH ENVIRONMENT.	

## EMC DECLARATION

---

1) This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS;

This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.

### BODY TEMPERATURE

ASTM laboratory accuracy requirements in the display range of 37° C (98° F to 102° F) for IR thermometers is  $\pm 0.2^{\circ}\text{C}$  ( $\pm 0.4^{\circ}\text{F}$ ), whereas for mercury in-glass and electronic thermometers, the requirement per ASTM Standards E667 -86 and E1112-86 is  $\pm 0.1^{\circ}\text{C}$  ( $\pm 0.2^{\circ}\text{F}$ ).

Body temperature varies from person to person and fluctuates during the course of the day. For this reason, it is suggested to know one's normal, healthy forehead temperature to correctly determine the temperature. Also it is for consumer use only, not for professional use.

- Body temperature runs approximately from 35.5° c to 37 .8° C (95.9° F ~ 100° F). To determine if one has a fever, compare the temperature detected with the person's normal temperature. A rise over the reference body temperature of 1° C (1 ° F) or more is generally indication of fever.

- Different measurement sites (rectal, axillary, oral, frontal, auricular) will give different readings. Therefore it is wrong to compare the measurement taken from different sites.

- Here below are typical temperatures for adults, based on different measurement sites:

- Rectal: 36.6° C to 38° C / 97.9° F ~ 100.4° F

- Axillary: 34.7° C to 37.3° C / 94.5° F ~ 99.1° F

2)\* Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.

3)\* Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!

4)\* Caution: this machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

5)\* The infrared body thermometer is suitable for hospitals except for near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.

6)\* Performance of the infrared body thermometer such as measure accuracy, information display would be lost or degraded due to EM DISTURBANCES.

7)\* Caution: the EMISSIONS characteristics of this equipment make it suitable for use in hospitals and a residential environment (CISPR 11 Class B).

### Guidance and manufacturer's statement- Electromagnetic emission

The Infrared body Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared body Thermometer should assure that it is used in such an environment.

Emission test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions IEC61000-3-3	Not applicable

### Guidance and manufacturer's declaration -Electromagnetic Immunity

The Infrared body Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared body Thermometer should assure that it is used in such an environment.

Anti-interference detection	IEC 60601-1 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 KV contact ±2 Kv, ±4 kV, ±8 kV, ±15 kV air	Same as the left
Electrical fast transient/burst IEC 61000-4-4	The input a.c. power ports: ±2 KV The input d.c. power ports: ±2 KV Signal input/output ports: ±1 KV	Not applicable
Surge IEC 61000-4-5	Input power ports: +0.5, +1.0 KV Signal input/output: +2.0 KV	Not applicable
Voltage interruption IEC 61000-4-11	0.5 cycles for > 95% (sync angle degrees):0, 45, 90, 135, 180, 225, 270, 315) 1 cycles for >95% UT (sync angle degrees):0) 250 (50Hz)/300 (60Hz) cycles for >95% UT (sync angle (degrees):0)	Not applicable
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m

NOTE UT is the a.c. mains voltage prior to application of the test level

## Guidance and manufacturer's declaration -Electromagnetic Immunity

Test Frequency (MHz)	Band (a) (MHz)	Service (a)	Modulation (b)	Maximum Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation (b) 18 Hz	1.8	0.3	27
450	430-470	GMRS460 FRS 460	FM (c) $\pm 5$ kHz deviation 1 kHz sine	2	0.3	28
710	704-787	LTE Band 13,17	Pulse modulation (b) 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation (b) 18 Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation (b) 217 Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation (b) 217 Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation (b) 217 Hz	0.2	0.3	9
5500						
5785						

(a) For some services, only uplink frequencies are included.

(b) The carrier shall be modulated using a 50% duty cycle square wave signal.

(c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

### Guidance and manufacturer's declaration -Electromagnetic Immunity

The Infrared body Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared body Thermometer should assure that it is used in such an environment.

Conducted RF  
IEC  
61000-4-6

3 Vrms  
150 kHz to 80 MHz

Not applicable

Radiated RF  
IEC  
6100-4-3

Professional healthcare  
environment: 3V/m  
Home healthcare environment:  
10 Vm  
80 MHz to 2700 MHz

Same as the left

NOTE 1: At 80 MHz and 800MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land-mobile radios, amateur radio, AM and Fm radio broadcast and TV broadcast cannot be predicated theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Infrared Body Thermometer is used exceeds the applicable RF compliance level above, the Infrared Body Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Infrared Body Thermometer.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### WARRANTY - TERMS AND CONDITIONS

This product benefits of a warranty of 24 months on material and manufacturing defects, starting from the date of purchase (see sales slip). The warranty of 24 months does not include damages caused by the usual use of parties identified as "consumable" (e.g., batteries, brush heads or parts subject to usury). The legal guarantee of 24 months is void if: 1. The product has undergone aesthetic damage due to improper use not in accordance with instructions in the manual. 2. This product has been modified and/or tampered with. 3. The cause of the failure was due to poor maintenance of the individual components and/or accessories and/or supplies (e.g. oxidation and/or scaling due to the retention of water or other liquids, sediment blocking the sensor, leak of corrosive liquid from batteries). The following is excluded from the legal guarantee of 24 months: 1. Costs related to replacements and/or repair of parts subject to wear or costs for ordinary maintenance of the product. 2. The costs and risks involved in transporting the product to and from the store where you purchased or otherwise authorized collection centre to receive the products under warranty. 3. Damage caused by or resulting from improper installation or improper use not in accordance with the directions in the instruction manual. 4. Damage due to natural disasters, accidental events or adverse conditions not compatible with the product. 5.

Defects that have a negligible effect on product performance. The manufacturer, distributor and all the parties involved in the sale do not assume any liability for losses and economic damage from any malfunction of the product. In accordance with current regulations the manufacturer, distributor and all the parties involved in the sale are not responding in any case for damages, including direct, indirect ones, loss of net income, loss of savings and additional damage and other details consequences going beyond the damage caused by the breach of warranty, contract, strict liability, wrongdoing or due to other causes, resulting from the use or inability to use the product and/or paper and electronic documents, including the lack of service. For further information on the help service visit the website **[www.nuvitababy.com](http://www.nuvitababy.com)**

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All rights reserved
Anteprima Brands International Ltd